

Marmo Canaloni s.r.l
Mr. Aldo Guadagni
54033 Carrara (MS) – Via Piace, 7A
ITALY

Testing of Natural stone

(4 enclosures)

Commission

Testing and assessment of marble to be used as outdoor façade panels. The marble is named *Canaloni*.

Samples and testing

Prior to the sampling, Björn Schouenborg, CBI and Bent Grek, Ramböll visited the quarries in Carrara 6-8th April for meetings with quarry owners and to obtain hand specimens for screening of potential marble types.

The subsequent sample preparation and transport to CBI was arranged by the quarry owner. CBI had provided detailed sampling and sample marking information but has no information about the actual sampling procedure. The requested number of specimens were received for discolouration test, laboratory bowing, petrography (AGA-Adjacent Grain Analysis) and flexural strength tests. More information about samples and tests is given in table 1.

Table 1. Overview of samples, tests and date of testing.

Sample-identity	No. of specimens	Tested	Property	Method	Test finalised
Canaloni	10	5	Laboratory bow test	EN 16306:2012	2016-08-09
	12	10	Flexural strength	EN 12372:2006	2016-06-15
	5	5	Flexural strength after bow test	EN 12372:2006	2016-08-09
	16	6	Discolouration test	EN 16140:2011	2016-07-18
	1	1	AGA Including petrographic analysis	EN 16306:2012 (Enclosure C) EN 12407:2007	2016-08-10

Test results

Average test results for the received specimens are given in the table below and detailed results in enclosures below.

Table 2. Test results.

Characteristics	Canaloni
Bowing potential (mm/m)	0,39
Flexural strength, reference (MPa), LEV *	14,5
Flexural strength, exposed (MPa), LEV *	11,2
Strength loss in % based on average values	26
Discoloration	T1
AGA	10

* Lower expected value (similar to characteristic value)

The Swedish Cement and Concrete Research Institute (CBI) Concrete and Stone

Performed by

Examined by

Lovise Sjöqvist

Björn Schouenborg

MSc Geology

PhD Mineralogy & Petrology

Enclosures

- 1 Determination of the resistance of marble to thermal and moisture cycling
- 2 Flexural strength of reference specimens and exposed specimens
- 3 Discolouration test
- 4 Petrographic description focussing on the AGA (Adjacent Grain Analysis) test: A quantitative measurement of the microstructure of the marble used as an indicative test of bowing potential and for follow up in quality control

Enclosure 1

Laboratory bow test for marble according to EN 16306:2012

The test was conducted by Lovise Sjöqvist, CBI, in the period between 13th July and 2nd August 2016.

Table 1. Average bowing of Canaloni marble.

	Average Bowing (mm/m)	Standard Deviation	Recommended limiting value (mm/m)	No. of Cycles
Canaloni	0,39	0,039	0,4	50

Figure 1. Bowing development during the 50 thermal and moisture cycles.





Figure 2. Bowing development during the 50 thermal and moisture cycles. Average values are shown. The recommended limit value is shown in green.



Enclosure 2

Flexural strength according to EN 12372 of reference specimens and exposed specimens

Reference samples

Method:	SS-EN 12372
Diameter of the loading device:	20 mm
Distance of the rollers:	350 mm
Load velocity:	0,25 MPa/s
Temperature:	40

Commission no:	6P00238
Sample:	Canaloni
Trading name:	Canaloni
Stone:	Marble
Date:	2016-06-15
Operator:	Lovise Sjöqvist

Specimen	Length, l (mm)	Width b (mm)	Height, h (mm)	Distance (mm)	Load (kN)	Flexural strength (MPa)
1	400	100.7	29.7	350	3.70	21.9
2	400	100.6	29.8	350	3.50	20.6
3	400	100.6	29.7	350	3.48	20.6
4	400	100.6	29.5	350	2.92	17.5
5	400	100.6	29.7	350	2.88	17.0
6	400	100.5	30.2	350	3.20	18.4
7	400	100.6	29.9	350	3.58	20.9
8	400	100.6	31.4	350	2.97	15.8
9	400	100.6	30.0	350	2.80	16.2
10	400	100.6	29.9	350	3.04	17.7
Mean value	400.0	100.6	30.0		3	18.7
Standard deviation	0.0	0.1	0.5		0.3	2.2
Variation coefficient					10.3	12
Characteristic value K	5%-fractile (75% confidens level)					14.5

Samples exposed to the “cladding bow test”

Natural Stone : Flexural strength 3-point load

Method:	SS-EN 12372
Diameter of the loading device:	20 mm
Distance of the rollers:	350 mm
Load velocity:	0,25 MPa/s
Temperature:	40

Commission no:	6P00238
Sample:	Canaloni
Trading name:	Canaloni
Stone:	Marble
Date:	2016-08-17
Operator:	Lovise Sjöqvist

Specimen	Length, l (mm)	Width b (mm)	Height, h (mm)	Distance (mm)	Load (kN)	Flexural strength (MPa)
1	400	100.5	29.6	350	2.35	14.0
2	400	100.6	29.7	350	2.37	14.0
3	400	100.5	29.7	350	2.46	14.6
4	400	100.7	29.7	350	2.14	12.7
5	400	100.6	29.7	350	2.56	15.1
6	400	100.8	29.6	350	1.94	11.5
7	400	100.5	30.1	350	2.28	13.2
8	400	100.6	30.0	350	2.60	15.1
Mean value	400.0	100.6	29.8		2	13.8
Standard deviation	0.0	0.1	0.2		0.2	1.2
Variation coefficient					9.4	9
Characteristic value K	5%-fractile (75% confidens level)					11.2

This means approximately 26 % reduction in flexural strength due to the bow test.

Enclosure 3

Discolouration test for marble according to EN 16140:2011, Annex B

The discolouration test has been found especially useful to assess the potential risk of brownish staining of marble. Given that there is any easily oxidized sulphides in the marble, they will be affected by repeated heating and cooling. The latter in a weak alkaline solution.

Figures a and b illustrate the test specimens before and after exposure to 20 cycles. **N.B!** The colours do not completely match the original, despite white balancing of the photos. The comparison to reference specimens (R1 and R2) are therefore crucial for the assessment!

Figure 1a. Test specimens before the exposure, including two reference specimens.





Figure 1b. Test specimens after the exposure, including two reference specimens that are not exposed.



Petrographic analysis with a focus on microstructure: Adjacent Grain Analysis (AGA) EN 12407:2007 and EN 16306:2012 Enclosure C

Test procedure

The microstructure is quantified using linear traverses for the grain size and the Adjacent Grain Analysis (AGA) for the complexity of the microstructure. Increased complexity is favourable as regards durability. This was done on thin section under a polarising microscope. The principle of the AGA-method is that the amount of adjacent grains is counted around the median-sized minerals/grains in the sample (Fig 1).




Figure 1. Schematic image of an igneous polygonal texture (left), and a seriate interlobate texture (right). The image also shows how the adjacent grains AG are counted.

Calcite crystal belongs to the hexagonal crystal system. In an “ideal” even-grained igneous texture, all calcite crystals share grain boundaries with six grains, which are referred to as *adjacent grains* AG (Figure 1). N.B! This is not ideal for the durability aspects! An increasingly irregular grain boundary or a more heterogeneous grain-size distribution will result in an increased number of adjacent grains. In a calcite marble with a heterogeneous grain size distribution, the largest crystals have the highest number of adjacent grains whereas the smallest crystals can have less than six adjacent grains. A larger number of AG occurs also in calcite marble with complex grain boundaries.

Enclosure 4

Test results

The analysed sample is a fine grained (0,02 – 0, 2mm) calcitic marble (Figure 1) with irregular grain boundaries and heterogeneous grain size (Figure 2), where a large number are calcite crystals, smaller than 0,1 mm due to recrystallization of the marble. For macroscopic images, see Enclosure 3. The median grain size is ~0.18 mm. The median of adjacent grains counted is 10 (median value) and 9.8 as a mean value (Figure 3).




Figure 1. The photograph shows the microstructure of the grains as well as grain size. Note smaller grains (~0,05 mm) along grain boundaries. Note also scale bar at the top left of the photo.

Enclosure 4




Figure 2. Grain-size distribution of the analysed sample.




Figure 3. The results of the adjacent grain analysis (AGA)

Enclosure 4

Comments to the results

According to results from bowing measurements on marble claddings in the laboratory and from facades, the marble types with a granoblastic microstructures have the highest bowing potentials. The normal range of AG in calcite marbles are between 6 and 10, where the marbles with 6 AG show the greatest bowing potential. This sample has **10 AG**. In comparison with earlier results, it is our experience that it should not have a potential to bow if used as a façade cladding, or it will be minimal.

CBI Swedish Cement and Concrete Research Institute

Materials, Borås

Magnus Döse

MSc. Geology

PRELIMINARY